



CLINICAL AND
LABORATORY
STANDARDS
INSTITUTE

8th Edition

CLSI PRE02™

Collection of Diagnostic Venous Blood Specimens

Sample

CLSI PRE02 provides procedures for the collection of diagnostic venous blood specimens, including line draws, blood culture collection, and venipuncture in children.

A standard for global application developed through the Clinical and Laboratory Standards Institute consensus process.

Collection of Diagnostic Venous Blood Specimens

Judith Dixon, MS, MT(ASCP), BS
Estelle Ninnemann, MLS(ASCP), MBA
Terri A. McElhattan, MHA, MLS, PBT(ASCP), CPI(ACA)
Michael O'Bryan, MD, MHA
Aurora Munoz Pedraza, MC EHDL

George F. Souza, PBT(ASCP), BS
Cherise Ens, MSc., MLT
Valerie Joyner, MBA, MLT(ASCP)
Charles E. Walker, MLT(ASCP)

Abstract

Clinical and Laboratory Standards Institute PRE02—*Collection of Diagnostic Venous Blood Specimens* provides a descriptive, stepwise process and procedures reflecting the quality system essentials format for diagnostic venous blood specimen collection. Special considerations for collections from vascular access devices, blood culture collection, and collections in isolation environments are included, as well as how to handle emergency situations. An expanded appendix section provides helpful tips for collecting specimens from pediatric and other challenging patients.

Clinical and Laboratory Standards Institute (CLSI). *Collection of Diagnostic Venous Blood Specimens*. 8th ed. CLSI standard PRE02. (ISBN 978-1-68440-270-0 [Print]; 978-1-68440-271-7 [Electronic]). Clinical and Laboratory Standards Institute, USA, 2025.

The Clinical and Laboratory Standards Institute consensus process, which is the mechanism for moving a document through two or more levels of review by the health care community, is an ongoing process. Users should expect revised editions of any given document. Because rapid changes in technology may affect the procedures, methods, and protocols in a standard or guideline, users should replace outdated editions with the current editions of CLSI documents. Current editions are listed in the CLSI catalog and posted on our website at www.clsi.org.

If you or your organization is not a member and would like to become one, or to request a copy of the catalog, contact us at:

P: +1.610.688.0100 **F:** +1.610.688.0700 **E:** customerservice@clsi.org **W:** www.clsi.org

Copyright ©2025 Clinical and Laboratory Standards Institute. Except as stated below, any reproduction of content from a CLSI copyrighted standard, guideline, or other product or material requires express written consent from CLSI. All rights reserved. Interested parties may send permission requests to permissions@clsi.org.

CLSI hereby grants permission to each individual member or purchaser to make a single reproduction of this publication for use in its laboratory procedures manual at a single site. To request permission to use this publication in any other manner, e-mail permissions@clsi.org.

To read CLSI's full Copyright Policy, please visit our website at <https://clsi.org/terms-of-use/>.

Suggested Citation

CLSI. *Collection of Diagnostic Venous Blood Specimens*. 8th ed. CLSI standard PRE02. Clinical and Laboratory Standards Institute; 2025.

Previous Editions:

August 1977, February 1979, March 1980, April 1984, July 1991, June 1998, December 2003, October 2007, April 2017

CLSI PRE02-Ed8

ISBN 978-1-68440-270-0 (Print)

ISBN 978-1-68440-271-7 (Electronic)

ISSN 1558-6502 (Print)

ISSN 2162-2914 (Electronic)

Volume 45, Number 3

Contents

Abstract	i
Committee Membership	iii
Foreword	vii
Chapter 1: Introduction	1
1.1 Scope	2
1.2 Standard Precautions	2
1.3 Terminology	2
Chapter 2: Path of Workflow	5
2.1 Patient Is Registered and Identified	7
2.2 Collection Test Request Is Received	7
2.3 Patient Is Approached and Greeted	7
2.4 Patient Is Identified	8
2.5 Precollection Requirements Are Assessed	9
2.6 Hands Are Cleansed	10
2.7 Patient Is Assessed	10
2.8 Patient Is Positioned	10
2.9 Specimen Is Collected	11
2.10 Specimen Is Labeled	22
2.11 Postvenipuncture Care Is Provided	23
2.12 Specimen Is Handled and Transported	24
Chapter 3: Blood Specimen That Cannot Be Obtained	25
Chapter 4: Complications	29
4.1 Accidental Arterial Puncture	30
4.2 Nerve Injury	30
4.3 Hematoma	30
4.4 Hemolysis	31
4.5 Monitoring Blood Volume Collected	31
4.6 First Aid	31
4.7 Dizziness, Syncope, or Unexpected Nonresponsiveness	32
4.8 Nausea	33
4.9 Vomiting	33
4.10 Convulsive Seizures	33
4.11 Incident Reports	33

Contents (Continued)

Chapter 5: Special Situations	35
5.1 Blood Culture Specimens	36
5.2 Therapeutic Drug Monitoring	37
5.3 Vascular Access Devices and Infusions	37
5.4 Patient Isolation	39
Chapter 6: Quality Management System Elements	41
6.1 Customer Focus	42
6.2 Facilities and Safety Management	42
6.3 Personnel Management	44
6.4 Supplier and Inventory Management	44
6.5 Documents and Records Management	47
6.6 Information Management	47
6.7 Nonconforming Event Management	47
6.8 Assessments	48
6.9 Continual Improvement	48
Chapter 7: Conclusion	49
Chapter 8: Supplemental Information	51
References	52
Appendix A. Difficult Collection	62
Appendix B. Sources for Establishing Maximum Blood Volumes to Be Collected From Patients Susceptible to Iatrogenic Anemia	65
Appendix C. Collections Proximal to an Intravenous Infusion Site	66
Appendix D. Preventing Syncope	67
The Quality Management System Approach	70

Foreword

During the collection and handling of blood specimens, numerous errors can occur that pose significant and avoidable risks to the patient and the health care professional. When global standards are not implemented, it is more likely that patients will be injured during the procedure, biologically representative specimens will not be obtained from patients, and test results will not be comparable from one facility to another or between visits for the patient.

The process and procedures detailed in CLSI PRE02 are intended to prevent specimen collection errors that threaten specimen quality, protect health care professionals from accidental exposure, and protect patients from the injuries, complications, and medical mistakes that can result from improperly collected specimens.

Since 1977, CLSI has recognized the importance of the preexamination phase of laboratory testing, including correct blood specimen collection and handling. Highly sophisticated testing technology cannot produce a good result from a poorly collected specimen.

Overview of Changes

CLSI PRE02-Ed8 replaces CLSI GP41-Ed7, published in 2017. Several changes were made in this edition, including:

- Removing duplicate information (ie, that included in CLSI PRE01¹)
- Harmonizing content with other CLSI preexamination documents
- Updating references

NOTE: The content of CLSI PRE02 is supported by the CLSI consensus process and does not necessarily reflect the views of any single individual or organization.

KEY WORDS

antecubital anatomy

blood specimen

complications from
phlebotomy

patient identification

phlebotomist

phlebotomy

veins

venipuncture

Chapter 1

Introduction

Sample

Collection of Diagnostic Venous Blood Specimens

1 Introduction

1.1 Scope

CLSI PRE02 establishes criteria for venous blood specimen collection for medical laboratory testing. These procedures are intended as an appropriate model for adoption by all health care providers responsible for blood specimen collection in outpatient and inpatient settings.

1.2 Standard Precautions

Because it is often impossible to know what isolates or specimens might be infectious, all patient and laboratory specimens are treated as infectious and handled according to “standard precautions.” Standard precautions are guidelines that combine the major features of “universal precautions and body substance isolation” practices. Standard precautions cover the transmission of all known infectious agents and thus are more comprehensive than universal precautions, which are intended to apply only to transmission of bloodborne pathogens. Published guidelines are available that discuss the daily operations of diagnostic medicine in humans and animals while encouraging a culture of safety in the laboratory.² For specific precautions for preventing the laboratory transmission of all known infectious agents from laboratory instruments and materials and for recommendations for the management of exposure to all known infectious diseases, refer to CLSI M29.³

1.3 Terminology

CLSI, as a global leader in standardization, is firmly committed to achieving global harmonization whenever possible. Harmonization is a process of recognizing, understanding, and explaining differences while taking steps to achieve worldwide uniformity. CLSI recognizes that medical conventions in the global metrological community have evolved differently in different countries and regions and that legally required use of terms, regional usage, and different consensus timelines are all important considerations in the harmonization process. CLSI recognizes its important role in these efforts, and its consensus process focuses on harmonization of terms to facilitate the global application of standards and guidelines. Table 1 is provided to clarify the intended interpretations of common terms.

Table 1. Common Terms or Phrases With Intended Interpretations

Term or Phrase	Intended Interpretation
“Needs to” or “must”	Explains an action directly related to fulfilling a regulatory and/or accreditation requirement or is indicative of a necessary step to ensure patient safety or proper fulfillment of a procedure
“Require”	Represents a statement that directly reflects a regulatory, accreditation, performance, product, or organizational requirement or a requirement or specification identified in an approved documentary standard
“Should”	Describes a recommendation provided in laboratory literature, a statement of good laboratory practice, or a suggestion for how to meet a requirement

2.8.1.2 Home Collections

Patients whose specimens are collected in their homes, where chairs with safety features are not available, must be placed in a recumbent position or in a chair with arm supports.

2.8.2 Recumbent Position

Health care professionals must position patients on their backs in a comfortable position by reclining the bed or phlebotomy chair or having patients lie flat after verification with the health care provider. For patients with whom communication is difficult (eg, sedated, cognitively impaired, language barrier), the health care professional should recline the patient with assistance, if necessary.


2.8.3 Pediatric Positioning

The pediatric patient can be securely positioned in the guardian's lap, or alternatively in a recumbent position on a table.

For infants, who are under 1 year of age and usually preambulatory, a heelstick (by capillary puncture) can be considered (refer to CLSI GP42¹² for information on capillary blood collection).

2.9 Specimen Is Collected

Figure 2 outlines the specimen collection process.



^a Four basic symbols are used in this process flow chart: oval (signifies the beginning or end of a process), arrow (connects process activities), box (designates process activities), diamond (includes a question with alternative “Yes” and “No” responses).

Figure 2. Specimen Collection Process Flow Chart^a

Sample



CLINICAL AND
LABORATORY
STANDARDS
INSTITUTE.

PRINT ISBN 978-1-68440-270-0

ELECTRONIC ISBN 978-1-68440-271-7

CLSI PRE02-Ed8